

MPR Bone Lever, U-Curved

MPR1058



Description

A bone lever with a distinctive “U”-curve design for soft tissue and bone retraction in hip, knee, and fracture repair surgeries. The shallow tip reduces risk of deep tissue injury, and the curved blade offers superior leverage and visibility in confined spaces.

Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization. Designed for optimal performance with minimal risk of trauma.

Made in USA

Product Specifications

Overall Length	7.28" (18.5cm)
Handle Length	4.33" (11cm)
Handle Width*	1.33" (3.38cm)
Tip Width	0.6" (1.6cm)

*measured at midpoint

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.