


## Stackable Step Stool

MPR140



### Description

Step stools for use in operating rooms. Designed to interlock via the tabs on the sides and stack to provide the exact height and platform size needed, and take up minimal storage space. Made from welded stainless steel to ensure strength and durability, with slip-resistant rubber tread surfaces and feet.

Made in USA.

### Product Specifications

Length	12" (30.48cm)
Width	18" (45.72cm)
Height	5" (12.7cm)

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).