



MPR Thin Glenoid Retractors



MPR1910



MPR1920

Description

Retractors used during shoulder arthroscopic surgeries to retract the humeral head and access the glenoid. Their thin construction compared to other retractors keeps them lightweight and maximizes precision and effectiveness, and the dual-pointed tip ensures a secure grip on tissue.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

Product Specifications

| # | DESCRIPTION | LENGTH | BLADE L. | BLADE W. | TIP L. | TIP W.* |
|---------|-------------|---------------|---------------|----------------|--------------|----------------|
| MPR1910 | Narrow | 11" (27.94cm) | 2.36" (6cm) | 0.51" (1.3cm) | 0.39" (1cm) | 0.48" (1.22cm) |
| MPR1920 | Wide | 10.24" (26cm) | 0.98" (2.5cm) | 0.88" (2.24cm) | 0.47" (1.2m) | 0.88" (2.24cm) |

*distance between points

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.