


MPR Reamer Hole Gauge

MPR210599000



Description

A measurement tool for checking and verifying reamer size, detect wear or undersizing on reamers, and ensure accurate canal sizing when fitting implants. The durable stainless steel construction minimizes wear on the gauge, ensuring consistently reliable measurements with years of repeated use, and is easy to clean.

Made in USA.

Product Specifications

Length	11" (27.94cm)
Width	4.75" (12.07cm)
Minimum Hole Ø	0.31" (8mm)
Maximum Hole Ø	0.89" (22.5mm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.