

## MPR Bone Spreader

MPR398.93



### Description

An instrument for spreading bone or bone fragments, providing better access and visibility during osteotomy and fracture repair procedures. The working end features one double-pronged forked tip and one single-pronged tip, allowing for insertion into tight spaces. The fluid action and speed-lock mechanism ensure controlled spreading of bone tissue, minimizing risk of trauma to the patient.

Made from solid stainless steel for durability and longevity, and easy cleaning/sterilization.

### Product Specifications

Overall Length	4.88" (12.38cm)
Handle Length	3" (7.6cm)
Handle Width	0.43" (1.1cm)
Tip Width	0.12" (3mm)
Beak Width	0.28" (7mm)
Maximum Spread	0.5" (12.5mm)
Curvature Angle	30°

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).