

## MPR Slap Hammer

MPR831



### Description

A slap hammer arm for use with MPR Orthopedics' curved jaw, straight jaw, and/or long nose vice grip locking pliers. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

NOT intended to work with MPR Flexible Osteotome Handles.

Made in USA.

### Product Specifications

Overall Length	14.75" (37.47cm)
Weight Length	3.75" (9.53cm)
Weight Diameter	1.25" (3.18cm)
Screw Length	0.5" (1.27cm)
Screw Diameter*	0.25" (6.35mm)
Overall Weight	1.9 lbs. (0.86kg)

\*including threads

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).