

## MPR Flexible Osteotome Handle, Bayonet Style

MPR948



### Description

A handle for MPR Flexible Osteotome Blades that use a bayonet style connector. Features a knurled grip, a faceted design to prevent rolling when placed down, and a quick-release button. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

Overall Length	6" (48.26cm)
Grip Length	3.5" (8.89cm)
Grip Diameter	0.64" (1.62cm)

### Single-Use Blades

#	LENGTH	WIDTH	END
MPR954	3" (7.62cm)	0.31" (8mm)	Straight
MPR955	3" (7.62cm)	0.39" (10mm)	Straight
MPR956	3" (7.62cm)	0.47" (12mm)	Straight
MPR957	3" (7.62cm)	0.79" (20mm)	Straight
MPR958	3" (7.62cm)	0.47" (12mm)	Curved
MPR960	5" (12.7cm)	0.31" (8mm)	Straight
MPR961	5" (12.7cm)	0.39" (10mm)	Straight
MPR959	5" (12.7cm)	0.79" (20mm)	Curved

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).