

## MPR Obturator

10089-1



### Description

An obturator with a ring handle for use with side-loading bone grafters. The obturator helps guide the grafter through soft tissue or into the joint capsule more easily than a blunt plunger. Made from stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Designed by Dr. Lanny Johnson, M.D.

### Product Specifications

Overall Length	8.5" (21.59cm)
Ring Diameter	1.25" (3.18cm)
Shaft Diameter	0.3" (0.75cm)

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).