

MPR Straight Rasp, 15° Bend

10705-15



Description

A rasp for restimulating the blood supply for complete meniscal repair in both open and arthroscopic procedures. The knurled grip allows for precise control, and its unique tooth design allows it to cut while pushing and pulling. Also used in conjunction with a suture passer in the Meniscus Mender II system. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Designed by Dr. Lanny Johnson, M.D.

Product Specifications

Overall Length	9" (22.86cm)
Length After Bend	2" (5.08cm)
Diameter	0.12" (3mm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.