



MPR Buck Femoral Cement Restrictor Inserter

MPR-11-2428



Description

Used during hip arthroplasty procedures to insert the cement restrictor. Designed to ensure proper and secure placement and fixation of the restrictor, and features one threaded end and one knurled end. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	13.25" (33.66cm)
Threaded End Ø	0.25" (0.64cm)
Knurled End Ø	0.75" (1.91cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.