



## MPR Dura-Flex Multi-Use Foam Pads



MPR-165-DF



MPR-931888-DF



MPR-931889-DF

### Description

High-durability foam pads for MPR Orthopedics leg positioning and prepping products (see *PRODUCT SPECIFICATIONS*). Their soft elastomeric foam construction offers both comfort and stability, and their polymeric covering is wear-resistant and moisture-impervious, allowing for repeated use.

All MPR Orthopedics foam products are latex-free. Made in USA.

### Product Specifications

#	COMPATIBILITY	LENGTH	DIAMETER
MPR-165-DF	MPR Single Loop Leg Prepper (MPR-164)	15" (38.1cm)	2" (5.08cm)
MPR-931888-DF	MPR Radioluscent Popliteal Leg Positioner (MPR-931888)	11" (28cm)	2.5" (6.35cm)
MPR-931889-DF	MPR Stainless Popliteal Leg Positioner (MPR-931889)	12.25" (31.1cm)	1.75" (4.45cm)

### Cleaning & Maintenance

Wipe down thoroughly with a surface disinfectant between uses. Discard when pad shows signs of wear and/or during routine maintenance.

### For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).