

MPR Moore T-Extractors



MPR-6839-2-060



MPR-6839-2-080

Description

Corkscrew-style extractors for removing the femoral head during hip arthroplasty procedures. The T-shaped handle is ergonomically designed for a secure, comfortable grip, and the threaded tip drills directly into the femoral head, allowing for removal without affecting nearby tissue. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

#	LENGTH	HANDLE WIDTH	TIP LENGTH	TIP WIDTH*
MPR-6839-2-060	6" (15.24cm)	3.57" (9.06cm)	1.64" (4.17cm)	0.3" (0.77cm)
MPR-6839-2-080	8" (20.32cm)	3.55" (9.02cm)	1.57" (3.98cm)	0.3" (0.77cm)

*measured at widest point

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.