

MPR T-Handle

MPR980



Description

T-Handle for use with trephines and driver shafts/slap hammers, featuring an ergonomic design to ensure a secure, comfortable grip and a Zimmer/Hall-style connector with a spring-loaded locking sleeve. Made from a wear- and heat-resistant plastic and solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

| | |
|-----------------|-----------------|
| Overall Length | 6.13" (15.56cm) |
| Handle Width | 5.51" (14cm) |
| Handle Diameter | 1.12" (2.85cm) |
| Shaft Diameter | 0.5" (1.27cm) |



For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.