

## MPR Solution Stands



MPR1000



MPR1001



MPR1002



MPR1003

### Description

Solution stands with included bowls. Stands and bowls are made from 16-gauge stainless steel for durability and reliability, with rubber bumpers around the basin to prevent noise when rolling, and 3" scuff-resistant casters designed for maneuverability and quiet movement.

Stands are available in either single-bowl or double-bowl designs, and with lower shelves or an H-brace construction.

Made in USA.

### Product Specifications

#	CONSTRUCTION	BOWLS	HEIGHT	WIDTH	DEPTH
MPR1000	H-Brace	Single	34" (86.36cm)	18" (45.72cm)	15" (38.1cm)
MPR1001	Shelf	Single	34" (86.36cm)	18" (45.72cm)	15" (38.1cm)
MPR1002	H-Brace	Double	34" (86.36cm)	33" (83.82cm)	15" (38.1cm)
MPR1003	Shelf	Double	34" (86.36cm)	33" (83.82cm)	15" (38.1cm)

### For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).