

Carry Cart

Description

A cart for transporting and storing stackable step stools such as the MPR Stackable Step Stool (MPR140). Made from heavy-gauge solid stainless steel for durability and reliability with a reinforced handle, and 3" double ball bearing swivel casters for easy maneuvering.

Made in USA.

Product Specifications

Length	21" (53.34cm)
Width	14" (35.56cm)
Height	39" (99.06cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR145

