

MPR Fukuda Retractors



MPR1504



MPR1505



MPR224510004

Description

Retractors used during shoulder orthopedic surgeries to posteriorly retract the humeral head and expose the anterior glenoid. Their distinctive design features a T-handle to ensure precise control, and a hole in the blade that securely grips onto the humeral head.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

Made in USA.

Product Specifications

#	DESCRIPTION	LENGTH	BLADE W.	T-HANDLE	HOLE WIDTH
MPR1504	MPR Fukuda-Type Retractor, Small	7.5" (19.05cm)	1.25" (3.18cm)	3.95" (10.03cm)	1" (2.54cm)
MPR1505	MPR Fukuda-Type Retractor, Large	7.5" (19.05cm)	1.5" (3.8cm)	3.95" (10.03cm)	1.25" (3.18cm)
MPR224510004	MPR Large Modified Fukuda Retractor	6.5" (16.51cm)	1.58" (4cm)	3.95" (10.03cm)	1.59" (4.05cm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.