

## MPR Wiltse Gelpi Retractor

MPR255612000



### Description

A self-retaining retractor with a ratchet-lock mechanism for spreading soft tissue to access a surgical site. The 70°-bent arms offer deep access and the outward-turned sharp tips ensure a secure grip on tissue.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

### Product Specifications

Overall Length	11" (27.94cm)
Arm Length	7" (17.78cm)
Tip Length	0.5" (1.27cm)
Tip Spread (closed)	1.5" (3.81cm)
Max Spread	6.88" (17.46cm)

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).