

MPR Curved Jaw Vise Grips/Pliers



Description

Pliers used for a number of applications, including revision, extraction, and explanting procedures. The serrated jaws ensure a firm grasp, the design is largely flat for easy cleaning while still being comfortable to hold, and the smooth gliding locking mechanism enables easy one-handed opening and closing. All models feature two inserts for a slap hammer. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

#	LENGTH	JAW LENGTH	BACK INSERT	SIDE INSERT
MPR820	5" (12.7cm)	0.88" (2.22cm)	Yes	Yes
MPR822	8" (20.32cm)	1.38" (3.49cm)	Yes	Yes
MPR823	10" (25.4cm)	1.5" (3.81cm)	Yes	Yes

Product Compatibility

#	DESCRIPTION
MPR831	MPR Slap Hammer

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.