

MPR Closed Loop Extractor

MPR946



Description

A loop extractor for removing a one-piece prosthesis stem during revision and extraction procedures. The loop design fits neatly around the the head of the prosthesis and the angled tip provides leverage for the extraction process. A slap hammer may be attached via the quick connector on the shaft. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	5.87" (14.92cm)
Overall Width	2.17" (5.52cm)
Tip Length	1.9" (4.84cm)
Tip Width	1.04" (2.63cm)
Tip Angle	30°
Wire Diameter	0.23" (0.59cm)

Product Compatibility

#	DESCRIPTION
MPR950	MPR Femoral Extractor Slap Hammer

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.