

## MPR Tendon Stripper

MPR-10760



### Description

An instrument primarily used for semitendinosus harvesting. Features a balanced handle for a secure grip and high level of control, and a conical cutting end for effective rotation. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Designed by Dr. Lanny Johnson, M.D.

### Product Specifications

Overall Length	17" (43.18cm)
Handle Length	5" (12.7cm)
Shaft Ø	0.12" (0.3cm)
Inside Cutting Ø	0.24" (0.6cm)

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).