



# MPR Buck Femoral Cement Restrictor Inserter

MPR-11-2428



## Description

Used during hip arthroplasty procedures for cement restrictor insertion. Designed to ensure proper and secure placement and fixation of the restrictor, and features one knurled end and one threaded end. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

## Product Specifications

Overall Length	13.25" (33.66cm)
Knurled End Ø	0.75" (1.91cm)
Threaded End Ø	0.25" (0.63cm)
Tip Ø	0.15" (0.39cm)

## For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).