

## MPR Graft Harvesting Retractor

MPR-1420



### Description

A retractor for use during ACL repair and other such procedures to harvest bone-tendon-bone (BTB) grafts. The forked tip hooks over the superior pole of the patella to securely retract the surrounding tissue, and the windowed shaft provides an optimal degree of access for harvesting the tendon.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

Overall Length	8.38" (21.27cm)
Handle Length	2.75" (6.99cm)
Window Length	3" (7.62cm)
Window Width (Ext.)*	1.18" (3cm)
Window Width (Int.)*	0.75" (1.91cm)
Tip Length	1" (2.54cm)
Tip Width	0.75" (1.91cm)

\*measured at widest point

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).