

## MPR Single Loop Leg Prepper

MPR-164



### Description

Also known as the “Shepard’s Hook”, a positioning device for preparing a patient for leg procedures, that easily attaches to an operating table via a single 5/8” Clark socket. Its simple, versatile design also allows it to be used to run cables and/or tubes over the operating table like a surgical picket fence. Made from solid stainless steel for durability, longevity, and easy cleaning.

Comes with a reusable foam pad (MPR-165).

Made in USA.

### Product Specifications

Height	33" (83.82cm)
Width	13" (33.02cm)
Weight Cap.*	50lbs (22.68kg)

\*approximate

### Accessories

#	DESCRIPTION
MPR-165	Reusable Pad w/ End Cap
MPR-165BLU	Single-Use Pad w/ End Cap
MPR-166BLU	Single-Use Pad w/o End Cap
MPR-165-DF	Dura-Flex Multi-Use Pad

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).