



Standard Tissue Graspers – Rongeurs



MPR-1790-01



MPR-1790-02

Description

Rongeurs with toothed tips for securely grasping and/or removing tissue during intraarticular procedures. Their ergonomic pistol-grip handle design ensures precise control and the low-profile shaft minimizes obstruction and offers a high degree of maneuverability within tight joint spaces.

MPR Orthopedics offers standard tissue graspers in two different sizes. Both are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization. Made in USA.

Product Specifications

#	DESCRIPTION	LENGTH	SHAFT	TIP WIDTH	TIP SPREAD
MPR-1790-01	5" Standard Tissue Grasper	8" (20.32cm)	5" (12.7cm)	0.21" (5.25mm)	0.55" (14mm)
MPR-1790-02	9" Standard Tissue Grasper	12" (30.48cm)	9" (22.86cm)	0.21" (5.25mm)	0.55" (14mm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.