



## MPR Triangle Knee Positioners



MPR-193P



MPR-293P

### Description

Triangular positioning devices for use during knee operations to position the leg without the need for a fracture table. Made from welded aluminum for reduced weight and easy cleaning between uses.

MPR Orthopedics offers triangle knee positioners in two different sizes, both of which are compatible with the MPR Triangle Pad (*sold separately, see ACCESSORIES*).

Made in USA.

### Product Specifications

#	DESCRIPTION	HEIGHT	WIDTH	DEPTH
MPR-193P	Large Triangle Knee Positioner	15" (38.1cm)	7" (17.78cm)	10" (25.4cm)
MPR-293P	Small Triangle Knee Positioner	12" (30.48cm)	7" (17.78cm)	7" (17.78cm)

### Accessories

#	DESCRIPTION	HEIGHT	WIDTH	THICKNESS
MPR-194	Single-Use Triangle Pad	10.75" (27.31cm)	7" (17.78cm)	1.38" (5.08cm)

### For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).