

## MPR Bone File

MPR202485004



### Description

A double-ended file for carving and shaping bone during revision and implant fitting procedures. Its two tips offer a total of eight filing/rasping faces with different levels of roughness, serration types, and surface shapes, allowing for a high degree of precision. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

Overall Length	8.5" (21.59cm)
Handle Length	3.5" (8.89cm)
Handle Width	0.39" (1cm)
Tip Length (each)	1.5" (10.8cm)
Tip Width	0.31" (0.8cm)

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).