

Meyerding Skid

MPR-3622



Description

A retractor for use during shoulder and hip surgeries to retract the humeral or femoral head. Features a T-handle for a secure, comfortable grip and precise control, and a scoop-shaped blade for a firm grip on the bone being retracted.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Product Specifications

Overall Length	13.38" (33.97cm)
Handle Width	3.13" (7.95cm)
Handle Diameter	0.31" (0.8cm)
Blade Length	4.5" (11.46cm)
Blade Width	1" (2.54cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.