



MPR Mallets, Phenolic Handles



MPR-399.42

MPR-399.43

Description

Universal orthopedic mallets. The handles are made from a durable, heat- and wear-resistant fabric-reinforced phenolic resin, and ergonomically designed to provide a secure, comfortable grip. The sides of the solid stainless steel heads are faceted to ensure stable placement.

MPR Orthopedics offers two models with different sizes and weights to suit your needs and/or preference. The handles can be safely steam-sterilized up to 290°F/143°C.

Made in USA.

Product Specifications

#	DESCRIPTION	LENGTH	HANDLE	HEAD DIAMETER	WEIGHT
MPR-399.42	MPR Mallet, 1.25 lb.	7.63" (19.37cm)	4.5" (11.4cm)	1.25" (3.18cm)	1.25 lb (570g)
MPR-399.43	MPR Mallet, 1.7 lb.	7.75" (19.69cm)	4.5" (11.4cm)	1.5" (3.81cm)	1.7 lb (771g)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.