



MPR Moreland Carbide Punch

MPR-6020-25



Description

A punch for use in cement removal and hip revision procedures. Features a tungsten carbide tip for optimal effectiveness and wear-resistance, and a faceted Radel®* handle to prevent rolling when placed down. Made from stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA. Original design by Dr. John R. Moreland, M.D.

Product Specifications

Overall Length	13" (33.02cm)
Handle Length	5" (12.7cm)
Handle Diameter	1.18" (2.99cm)
Shaft Length	7.42" (18.85cm)
Shaft Diameter	0.23" (0.58cm)
Tip Length	0.45" (1.14cm)
Tip Diameter	0.15" (0.39cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

*Compared with other materials commonly used on handles for surgical instruments, Radel® is engineered for harsh environments, and can be cleaned, disinfected, and steam-sterilized over 1,000 cycles without significant loss of properties.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com