



## MPR Moore T-Extractors



MPR-6839-2-060



MPR-6839-2-080

### Description

Corkscrew-style extractors used to remove the femoral head during Total Hip Arthroplasty procedures. The T-shaped handle is ergonomically designed for a secure, comfortable grip, and the threaded tip drills directly into the femoral head, allowing for removal without affecting nearby tissue. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

| #              | LENGTH       | HANDLE WIDTH   | TIP LENGTH     | TIP WIDTH*    |
|----------------|--------------|----------------|----------------|---------------|
| MPR-6839-2-060 | 6" (15.24cm) | 3.57" (9.06cm) | 1.64" (4.17cm) | 0.3" (0.77cm) |
| MPR-6839-2-080 | 8" (20.32cm) | 3.55" (9.02cm) | 1.57" (3.98cm) | 0.3" (0.77cm) |

\*measured at widest point

### For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).