

MPR PCL Tibial Retractor

MPR-86-9050



Description

Also known as a “Pickle Fork” retractor, a retractor for use during knee procedures to retract the tibia away from the femur. Its two-pronged tip allows it to fit into the femoral condylar notch while straddling the posterior cruciate ligament, and the handle features three holes for attaching weights.

All MPR Orthopedics retractors are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and are designed for optimal performance with minimal risk of trauma.

Made in USA.

Product Specifications

Overall Length	11” (4.99cm)
Tip Length	1.11” (2.81cm)
Tip Width	0.85” (2.15cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.