



MPR Popliteal Leg Positioners



MPR-931888



MPR-931889



MPR-931890

Description

Knee supports for TKA procedures. All can attach to a standard OR table rail using the MPR Over-the-Rail Universal Clamp (MPR-ORUC), and can be draped during surgery, thus not requiring sterilization. A variety of single-use and limited-reuse pads are also available, minimizing risk of cross-contamination. Made in USA.

Product Specifications

#	DESCRIPTION	HEIGHT	BAR LENGTH	BAR Ø
MPR-931888	MPR Radiolucent Popliteal Leg Positioner	16.5" (41.91cm)	12" (30.48cm)	2" (5.08cm)
MPR-931889	MPR Stainless Steel Popliteal Leg Positioner	11.5" (29.21cm)	11.5" (29.21cm)	0.63" (1.59cm)
MPR-931890	MPR Popliteal Leg Positioner, Slim	16.5" (41.91cm)	11" (27.94cm)	2" (5.08cm)

Accessories

#	DESCRIPTION	COMPATIBILITY
MPR-931888BLU*	MPR Single-Use Radiolucent Leg Positioner Pad, Blue	MPR-931888, MPR-931890
MPR-TKSP-R†	MPR Black Conductive Single Use Pad	MPR-931888, MPR-931890
MPR-931888-DF	MPR Dura-Flex Multi-Use Radiolucent Leg Positioner Pad	MPR-931888, MPR-931890
MPR-931889BLU†	MPR Single Use Pad Blue	MPR-931889
MPR-TKSP-SS*	MPR Black Conductive Single Use Pad	MPR-931889
MPR-931889-DF	MPR Dura-Flex Multi-Use Pad	MPR-931889

*included

†pictured

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.