

## MPR Bone Pick/Awl

MPR-BP



### Description

The original orthopedic instrument: a tool for creating and/or enlarging holes in bone. MPR Orthopedics' take on this Paleolithic classic features an extended length compared to market standard designs, offering a greater degree of access. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

Overall Length	9.75" (24.77cm)
Handle Length	4.25" (10.8cm)
Shaft Diameter	0.2" (4.7mm)

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).