



MPR Flexible Osteotomes, Thin



Description

Thin flexible blades used in conjunction with a compatible handle for osteotomies and implant loosening. All are made from solid stainless steel for maximal effectiveness, and feature a circular-style connector.

Made in USA. Available either sterile or nonsterile. Intended for single use *only*; safely dispose after use.

Product Specifications

#	BLADE LENGTH	BLADE WIDTH
MPR575	1.5" (3.81cm)	0.24" (0.6cm)
MPR576	1.5" (3.81cm)	0.39" (1cm)
MPR574	3" (7.62cm)	0.24" (0.6cm)
MPR577	3" (7.62cm)	0.31" (0.8cm)
MPR578	3" (7.62cm)	0.39" (1cm)

#	BLADE LENGTH	BLADE WIDTH
MPR579	3" (7.62cm)	0.47" (1.2cm)
MPR580	3" (7.62cm)	0.79" (2cm)
MPR581	5" (12.7cm)	0.31" (0.8cm)
MPR582	5" (12.7cm)	0.39" (1cm)

Product Compatibility

#	DESCRIPTION
MPR599	MPR Osteotome Handle

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.