

## MPR K-Wire Cases



MPR-2071-10



MPR-424-80



MPR-424-85



MPR-424-90

### Description

Folding cases/racks for holding Kirschner Wires, commonly known as “K-Wires”, and Steinmann Pins. Improved from an older design, these cases offer up to more than triple the carrying capacity, while still being optimized for compactness and convenience, and are labeled with both US standard and metric units.

MPR Orthopedics offers cases in three different sizes, as well as a miniature rack with a convenient tilt-stand design and folding handle. All are compatible with all manufacturers’ K-Wires and Steinmann Pins, and are made from solid stainless steel to ensure durability and longevity.

Made in USA.

### Product Specifications

#	DESCRIPTION	LENGTH	WIDTH	HEIGHT*
MPR-2071-10	Mini K-Wire Rack	---	4.33" (11cm)	4.75" (12.07cm) <sup>†</sup>
MPR-424-80	Mini K-Wire Case	4.5" (11.43cm)	3.63" (9.22cm)	2.5" (6.35cm)
MPR-424-85	K-Wire Case	6.5" (16.5cm)	3.63" (9.22cm)	4.13" (10.49cm)
MPR-424-90	Large Pin/Wire Case	9.63" (24.46cm)	6.19" (15.72cm)	5.5" (13.97cm)

\*open/unfolded

<sup>†</sup>when handle is folded down

See reverse for information on carrying capacities.

### For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).

## Carrying Capacities

### MPR Mini K-Wire Rack

MPR-2071-10

The MPR Mini K-Wire Rack is designed to hold 4" (10.16cm) wires in the following quantities:

DIAMETER	#
0.028" (0.7mm)	42
0.035" (0.9mm)	39
0.045" (1.1mm)	48
0.054" (1.4mm)	24
0.062" (1.6mm)	22
<b>TOTAL:</b>	<b>175</b>

### MPR Mini K-Wire Case

MPR-424-80

The MPR Mini K-Wire Case is designed to hold 3" (7.62cm) and 4" (10.16cm) wires in the following quantities:

DIAMETER	#
0.028" (0.7mm)	19
0.035" (0.9mm)	21
0.045" (1.1mm)	25
0.054" (1.4mm)	15
0.062" (1.6mm)	25
<b>TOTAL:</b>	<b>105</b>

### MPR K-Wire Case

MPR-424-85

The MPR K-Wire Case is designed to hold 5" (12.7cm) and 6" (15.24cm) wires/pins in the following quantities:

DIAMETER	#
0.028" (0.7mm)	30
0.035" (0.9mm)	28
0.045" (1.1mm)	36
0.054" (1.4mm)	30
0.062" (1.6mm)	31
0.078" (2.0mm)	29
<b>TOTAL:</b>	<b>184</b>

### MPR Large K-Wire Case

MPR-424-90

The MPR Large K-Wire Case is designed to hold 9" (22.86cm) wires/pins in the following quantities:

DIAMETER	#
0.035" (0.9mm)	8
0.045" (1.1mm)	9
0.062" (1.6mm)	8
0.078" (2.0mm)	8
0.094" (2.4mm)	9
0.110" (2.8mm)	9
0.125" (3.2mm)	8
0.142" (3.6mm)	8
0.157" (4.0mm)	8
0.189" (4.8mm)	8
<b>TOTAL:</b>	<b>83</b>