



MPR Tissue Protection Sleeve

Description

Also known as a Freehand Tissue Protection Sleeve. An instrument comprising a handle and a barrel, used to protect soft tissue and as an instrument guide during fine reaming and drilling procedures. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	8" (20.32cm)
Handle Length	4.88" (12.38cm)
Barrel Length	3" (7.62cm)
Outer Diameter	0.41" (1.04cm)
Inner Diameter	0.32" (0.81cm)

MPR-TPRO



For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.