

MPR Capsule Retractor, 1 Prong

MPR1009



Description

A retractor for use during shoulder procedures to retract the glenoid capsule. The minimal, bent design provides a high degree of leverage, and the long, slender single-pronged tip enhances precision during use.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	10" (25.4cm)
Handle Length	3.5" (8.89cm)
Handle Bend	40°
Tip Length	1" (2.54cm)
Tip Angle	30°

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.