



Stainless Steel Linen Hampers



Description

Carts for transporting linens to and from the OR. Made from heavy gauge stainless steel for durability and reliability, with 3" double ball bearing swivel casters for easy maneuvering in operating rooms and hallways. Available in circular-top or triangle-top models.

Made in USA.

Product Specifications

CIRCULAR-TOP MODELS

#	HEIGHT	DIAMETER	CIRCUMFERENCE	BAG DIAMETER
MPR2000	35" (88.9cm)	25" (63.5cm)	78.54" (199.49cm)	25" (63.5cm)
MPR2001	35" (88.9cm)	18" (45.72cm)	56.55" (143.64cm)	18" (45.72cm)

TRIANGULAR-TOP MODELS

#	HEIGHT	WIDTH	PERIMETER	BAG DIAMETER
MPR2005	35" (88.9cm)	22.7" (57.66cm)	78.63" (199.72cm)	25" (63.5cm)
MPR2006	35" (88.9cm)	16.53" (41.99cm)	57.25" (145.42cm)	18" (45.72cm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.