

MPR Forked Retractor

MPR230786002



Description

A retractor for use during shoulder arthroplasty procedures for retracting the deltoid muscle and other soft tissues, featuring a fenestrated blade, a two-pronged tip to ensure a secure grip on tissue, and a hole in the handle for attaching weights.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Product Specifications

Overall Length	10.83" (27.5cm)
Handle Length	7.43" (18.86cm)
Blade Length	1.65" (4.2cm)
Blade Width	0.94" (2.39cm)
Tip Length	1.31" (3.33cm)
Tip Width*	0.72" (1.83cm)
Prong Width	0.8" (4.4mm)
Point Spread	0.56" (14.25mm)

*measured at widest point

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.