



MPR Blount Knee Retractor

MPR248500000



Description

A retractor used during knee procedures to retract the surrounding tissue such as tendons and ligaments. The broad leaf-shaped blade and long, narrow tip provide deep access and precision while enabling safe retraction of larger amounts of tissue, and the deeply curved, hook-shaped design maximizes leverage.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Product Specifications

Overall Length	7" (17.75cm)
Handle Length	4.94" (12.55cm)
Handle Width	0.59" (1.5cm)
Blade Length	1.02" (2.6cm)
Blade Width	0.69" (1.75cm)
Tip Length	1.18" (3cm)
Tip Width	0.3" (7.5mm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.