

## MPR MI Sciatic Nerve Retractor

MPR259807150



### Description

A retractor for use during hip procedures utilizing the minimally invasive technique to retract the sciatic nerve. Features a low-profile 90° bent design, a serrated round tip to ensure a secure grip on tissue, and a hole in the handle for attaching weights.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

Overall Length	7.99" (20.3cm)
Handle Length	5.14" (13.06cm)
Blade Length	1.65" (4.2cm)
Blade Width	1" (2.54cm)
Tip Length	0.65" (1.65cm)
Tooth Length	0.09" (2.34mm)
Tooth Width	0.1" (2.65mm)
Point Width	0.02" (0.5mm)

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).