

MPR Microplasty Ligamentum Cutter

MPR31-5555-37



Description

A cutter used to sever the ligamentum teres femoris during minimally invasive THA procedures, with a knurled handle to ensure a secure, comfortable grip, a thin, curved blade allowing access to the ligament without requiring the femoral head to first be dislocated from the acetabulum, and a specialized tip design that simultaneously tensions and cuts the ligament.

Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	10.5" (26.67cm)
Handle Length	5.5" (13.97cm)
Handle Diameter	0.75" (1.91cm)
Blade Length	5" (12.7cm)
Blade Width	0.75" (1.91cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.