



MPR Reusable Stainless Steel Finger Traps



MPR3280



MPR3285



MPR3290

Description

Finger traps for upper extremity traction made from stainless steel wire, allowing them to be cleaned and sterilized between uses to help cut down on single-use product waste. It is recommended that the patient's fingers are wrapped in sterile drapes before application to reduce skin trauma.

MPR Orthopedics offers three different sizes. All are made from solid stainless steel, ensuring durability and longevity, and allowing for easy cleaning/sterilization.

Made in USA.

Product Specifications

#	SIZE	LENGTH	MAX OPENING DIAMETER
MPR3280	Small	4.25" (10.8cm)	0.44" (1.11cm)
MPR3285	Medium	5.75" (14.61cm)	0.75" (1.91cm)
MPR3290	Large	5.25" (13.34cm)	1" (2.54cm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.