

## MPR Townley Femur Caliper

MPR3895



### Description

A caliper primarily used for measuring the femoral head. Features measurements in U.S. customary and metric units up to 4"/100mm and a thumbscrew for quickly loosening/tightening the slide. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

|                |                 |
|----------------|-----------------|
| Overall Length | 4.56" (11.59cm) |
| Handle Width   | 0.79" (2cm)     |
| Max Depth      | 1.39" (3.49cm)  |

\*measured at widest point

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).