



## MPR Picket Fence Leg Preppers



### Description

Limb positioning devices for use in patient prepping for leg surgery. These preppers attach easily to the operating table via either a single  $\frac{5}{8}$ " Clark socket, or a pair of blade-style sockets (see *PRODUCT SPECIFICATIONS*), and can be freely raised and/or lowered to achieve the required height.

MPR Orthopedics offers two different picket fence preppers, for either one or both legs. Both are made from solid stainless steel, with a double urethane coating, for durability, longevity, and easy cleaning/disinfecting. Made in USA.

### Product Specifications

#	DESCRIPTION	HEIGHT	OPERATIVE HEIGHT RANGE	WIDTH	ATTACHMENT
MPR611	Single	30" (76.2cm)	4"–27" (10.16cm–68.58cm)	19" (48.26cm)	Clark Socket
MPR611B	Double	33" (83.82cm)	4"–30" (10.16cm–76.2cm)	36" (91.44cm)	2 Blade Sockets

### Cleaning & Maintenance

Between uses, wipe positioner down with a hospital- and/or EPA-approved germicide, or a **10:1** bleach solution or quaternary ammonia solution, then rinse and let dry. Avoid using phenolic cleaners or abrasive materials.

### For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).