



## MPR Straight Jaw Vise Grips/Pliers



### Description

Pliers used for a number of applications, including revision, extraction, and explanting procedures. Their flat design is optimized for comfortable one-handed operation with a minimal profile, and features serrated jaws for a secure grip, and up to two inserts for attaching a slap hammer. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

#	LENGTH	JAW LENGTH	BACK INSERT	SIDE INSERT
MPR819	10" (25.4cm)	1.38" (3.49cm)	No	No
MPR824	7" (17.78cm)	1.38" (3.49cm)	Yes	Yes
MPR825	10" (25.4cm)	1.38" (3.49cm)	Yes	Yes

### Product Compatibility

#	DESCRIPTION
MPR831	MPR Slap Hammer

### For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).