


MPR Long Nose Vise Grips/Pliers



Description

Pliers used for a number of applications, including revision, extraction, and explanting procedures. Their flat design is optimized for comfortable one-handed operation with a minimal profile, and features long tapered jaws for precision and a secure grip, and up to two inserts for attaching a slap hammer. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

#	LENGTH	JAW LENGTH	BACK INSERT	SIDE INSERT
MPR826	6" (15.24cm)	2.5" (5.72cm)	Yes	No
MPR827	9" (22.86)	2.75" (6.99cm)	Yes	Yes
MPR828	11" (27.94cm)	2.5" (5.72cm)	Yes	Yes

Product Compatibility

#	DESCRIPTION
MPR831	MPR Slap Hammer

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.