


## MPR Crochet Hook



### Description

A hook used during revision procedures for cement removal/extraction, particularly removal of cement restrictors from the femur. Features an ergonomically designed Radel®\* handle to ensure a secure, comfortable grip, a long, thin shaft with a tapered end for deep access and precision, and a sharp, angular hooked tip. Made from stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

Overall Length	13" (33.02cm)
Handle Width	5.39" (13.7cm)
Shaft Length	3.75" (9.52cm)
Shaft Diameter	0.25" (6.4mm)
Tip Length	0.5" (1.27cm)
Tip Width	0.25" (6.4mm)

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).

\*Compared with other materials commonly used on handles for surgical instruments, Radel® is engineered for harsh environments, and can be cleaned, disinfected, and steam-sterilized over 1,000 cycles without significant loss of properties.