



MPR Rounded Acetabular Cement Splitter

MPR925



Description

A splitter used during hip extraction and cement removal/revision procedures to extract the cement mantle after a cemented acetabular component has been removed. The tip is rounded to closely match the curve of the extracted acetabular component, and the sharp edges allow for easy fragmentation and removal of the cement. Features an ergonomic Radel®* handle for a secure, comfortable grip. Made from stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	14" (35.56cm)
Handle Length	5.13" (13cm)
Shaft Length	8.63" (21.91cm)
Shaft Diameter	0.25" (6.4mm)
Tip Length	1.13" (2.86cm)
Tip Width	0.41" (1.03cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

*Compared with other materials commonly used on handles for surgical instruments, Radel® is engineered for harsh environments, and can be cleaned, disinfected, and steam-sterilized over 1,000 cycles without significant loss of properties.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com