

MPR Hook Stem Extractor

MPR947



Description

A hook-shaped extractor for extracting implants that feature a transverse extraction hole. The hook inserts easily into the hole, at which point a slap hammer is attached using the quick connector on the extractor to provide the extraction force. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	4.92" (12.5cm)
Overall Width	1.57" (4cm)
Hook Length	2.96" (7.53cm)
Tip Diameter	0.2" (0.52cm)
Tip Angle	65°

Product Compatibility

#	DESCRIPTION
MPR950	MPR Femoral Extractor Slap Hammer

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.